



x-tracta™ Disposable Gel Extraction Tool


- Eliminates cross contamination
- Minimizes UV exposure and damage
- Cuts uniform gel slices
- Replaces hazardous blades

The x-tracta agarose gel tool gives you cleaner DNA in less time and with less effort. This disposable device cleanly cuts and holds uniform gel slices, eliminating cross-contamination from razor blades, tweezers, or other tools. The process is quick, reducing damaging UV exposure. Gel bands can be released directly into microcentrifuge tubes for the next extraction step.

5454-2500 Pack of 25
5454-0100 Pack of 100



INSTRUCTIONS FOR USE



In accordance with good laboratory practice, use the x-tracta tool once, then discard it to prevent cross-contamination.
Note: Do not squeeze the bulb until step 4. At step 4, use a quick snap-action squeeze to release the gel.