

ErgoOne®

Test Specifications

| | Test Volume | Systematic Error (Inaccuracy) | Random Error (Imprecision) |
|-----------------------|--------------|-------------------------------|----------------------------|
| Single Channel | | | |
| 0.1 - 2.5 μ l | 0.1 μ l | \pm 24% | \pm 12% |
| | 1.25 μ l | \pm 2.5% | \pm 1.5% |
| | 2.5 μ l | \pm 1.4% | \pm 0.7% |
| 0.5 - 10 μ l | 1.0 μ l | \pm 2.5% | \pm 1.8% |
| | 5 μ l | \pm 1.5% | \pm 0.8% |
| | 10 μ l | \pm 1.0% | \pm 0.4% |
| 2 - 20 μ l | 2 μ l | \pm 5.0% | \pm 1.5% |
| | 10 μ l | \pm 1.2% | \pm 0.6% |
| | 20 μ l | \pm 1.0% | \pm 0.3% |
| 10 - 100 μ l | 10 μ l | \pm 3.0% | \pm 1.0% |
| | 50 μ l | \pm 1.0% | \pm 0.3% |
| | 100 μ l | \pm 0.8% | \pm 0.2% |
| 20 - 200 μ l | 20 μ l | \pm 2.5% | \pm 0.7% |
| | 100 μ l | \pm 1.0% | \pm 0.3% |
| | 200 μ l | \pm 0.6% | \pm 0.2% |
| 100 - 1000 μ l | 100 μ l | \pm 3.0% | \pm 0.6% |
| | 500 μ l | \pm 1.0% | \pm 0.2% |
| | 1000 μ l | \pm 0.6% | \pm 0.2% |
| 500 - 5000 μ l | 500 μ l | \pm 2.4% | \pm 0.6% |
| | 2500 μ l | \pm 1.2% | \pm 0.25% |
| | 5000 μ l | \pm 0.6% | \pm 0.15% |
| Multichannel | | | |
| 0.5 - 10 μ l | 1.0 μ l | \pm 8.0% | \pm 5.0% |
| | 5 μ l | \pm 4.0% | \pm 2.0% |
| | 10 μ l | \pm 2.0% | \pm 1.0% |
| 10 - 100 μ l | 10 μ l | \pm 3.0% | \pm 2.0% |
| | 50 μ l | \pm 1.0% | \pm 0.8% |
| | 100 μ l | \pm 0.8% | \pm 0.3% |
| 30 - 300 μ l | 30 μ l | \pm 3.0% | \pm 1.0% |
| | 150 μ l | \pm 1.0% | \pm 0.5% |
| | 300 μ l | \pm 0.6% | \pm 0.3% |

