


TipOne Repeating Tip Selection Guide for USA Scientific Dispenser

Volume Selection Dial					Tip Choice
1	2	3	4	5	
Sample Volume					
2 μ l	4 μ l	6 μ l	8 μ l	10 μ l	0.1 ml
10 μ l	20 μ l	30 μ l	40 μ l	50 μ l	0.5 ml
20 μ l	40 μ l	60 μ l	80 μ l	100 μ l	1.0 ml
50 μ l	100 μ l	150 μ l	200 μ l	250 μ l	2.5 ml
100 μ l	200 μ l	300 μ l	400 μ l	500 μ l	5.0 ml
200 μ l	400 μ l	600 μ l	800 μ l	1000 μ l	10.0 ml
500 μ l	1.0 ml	1.5 ml	2.0 ml	2.5 ml	25 ml
1.0 ml	2.0 ml	3.0 ml	4.0 ml	5.0 ml	50 ml
Number of Dispensing Steps					 www.usascientific.com
48	23	15	11	8	